

## Human IL-17RB protein (Recombinant) (C-His) (STJP020303)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

### TARGET INFORMATION

Gene ID	<a href="#">55540</a>
Gene Symbol	<a href="#">IL17RB</a>
UniProt ID	<a href="#">I17RB_HUMAN</a>
Immunogen Region	Arg18-Per292

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---