

Human KIR2DL5/CD158f1 protein (Recombinant) (C-His&Avi) (STJP020308)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	57292
Gene Symbol	KIR2DL5A
UniProt ID	K12LA_HUMAN
Immunogen Region	His22-His240

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---