

## Human LILRA4/ILT-7/CD85g protein (Recombinant) (C-His&Avi) (STJP020312)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

### TARGET INFORMATION

Gene ID	<a href="#">23547</a>
Gene Symbol	<a href="#">LILRA4</a>
UniProt ID	<a href="#">LIRA4_HUMAN</a>
Immunogen Region	Glu24-Asn446

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---