

## Cynomolgus IL-31 protein (Recombinant) (N-His) (STJP020323)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |
|---------------------|------------------------------------------------------------------------------------------------|

### TARGET INFORMATION

|                  |              |
|------------------|--------------|
| Immunogen Region | Leu27-Thr163 |
|------------------|--------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|