

## Human DHLAG/CD74 protein (Recombinant) (N-His) (STJP020332)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	------------------------------------------------------------------------------------------------

### TARGET INFORMATION

Gene ID	<a href="#">972</a>
Gene Symbol	<a href="#">CD74</a>
UniProt ID	<a href="#">HG2A_HUMAN</a>
Immunogen Region	Gln73-Met232

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------