

Human RNF43 protein (Recombinant) (C-hFC) (STJP020349)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54894
Gene Symbol	RNF43
UniProt ID	RNF43_HUMAN
Immunogen Region	Gly24-Tyr197
Immunogen Sequence	GFGRGTGLVLAIVESERSAE QKAIIRVIPLKMDPTGKLNLTLEGVFAGVAEITPAEGKLM QSHPLYLCNASDDDNLEPGF ISIVKLESPRRAPRPCLSLA SKARMAGERGASAVLFDITE DRAAAEQQLQQPLGLTWPVVL IWGNDAEKLMEFVYKNQKAHVRIELKEPPAWPDY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---