

Mouse CLEC9A/CD370 protein (Recombinant) (C-hFc) (STJP020351)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, pH 7.4.
Storage Instruction	Store at -70 °C. This product is stable at ≤ -70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	232414
Gene Symbol	Clec9a
UniProt ID	CLC9A_MOUSE
Immunogen Region	Lys57-Ile264
Immunogen Sequence	KFFQVSSLVLEQQRLIQQD TALVNLTQWQRKYTLEYCQA LLQRSLHSGTDASTGPIVLT SPQMVPTLDSKETGSDCSP CPHNWIQNGKSCYYVFERWE MWNISKKSCLKEGASLFQID SKEEMEFISSIGLKGNGKY WVGVFQDGISGSWFWDGSS PLSDLLPAERQRSAGQICGY LKDSTLISDKCDSWKYFICE KKAFGSCI
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.