

Human Thrombopoietin R/MPL/CD110 protein (Recombinant) (C-His) (STJP020353)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4352
Gene Symbol	MPL
UniProt ID	TPOR_HUMAN
Immunogen Region	Gln26-Trp491
Immunogen Sequence	QDVSLASDSEPLKCFSTRF EDLTCFWDEEEAAPSGTYQL LYAYPREKPRACPLSSQSMP HFGTRYVCQFPDQEEVRLFF PLHLWVKNVFLNQTRTQRVL FVDSVGLPAPPSIIKAMGGS QPGELQISWEEPAPFISDFL RYELRYGPRDPKNSTGPTVI QLIATETCCPALQRPHSASA LDQSPCAQPTMPWQDGPQT SPSREASALTAEGGSLISG LQPGNSYWLQLRSEPDGIS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---