

## Human Alkaline phosphatase (Placental type) /ALPP protein (Recombinant) (C-6His) (STJP020368)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris-HCl, 150mM NaCl, 0.2mM CaCl <sub>2</sub> , pH 7.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">250</a>
Gene Symbol	<a href="#">ALPP</a>
UniProt ID	<a href="#">PPB1_HUMAN</a>
Immunogen Sequence	IIPVEEENPDFWNREAAEAL GAAKKLQPAQTAAKNLIIFL GDGMGVSTVTAARILKGQKK DKLGPEIPLAMDRFPYVALS KTYNVDKHVPDSGATATAYL CGVKGNFQTIGLSAAARFNQ CNTTRGNEVISVMNRKAKAG KSVGVVTTTRVQHASPAGTY AHTVNRNWYSDADVPASARQ EGCQDIATQLISNMDIDVIL GGGGRKYMFRMGTPDPEYPDD YSQGGTRLDGKNLVQEWLA

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---