

Mouse S100-A9 protein (Recombinant) (N-His) (STJP020369)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as 0.22 µg/ml filtered solution in PBS, 1mM DTT, 20% Glycerol pH 7.5.
Storage Instruction	Store at -70 °C. This product is stable at ≤ -70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20202
Gene Symbol	S100a9
UniProt ID	S10A9_MOUSE
Immunogen Region	Ala2-Lys113
Immunogen Sequence	ANKAPSQMERSITTIIDTFH QYSRKEGHPDTLSKKEFRQM VEAQLATFMKKEKRNEALIN DIMEDLDTNQDNQLSFEECM MLMAKLIFACHEKLHENNPR GHGSHGKGCGK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---