

Mouse Galectin-3/LGALS3 protein (Recombinant) (C-His) (STJP020370)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Immunogen Region	Ala2-Ile264
Immunogen Sequence	ADSFSLNDALAGSGNPNPQG YPGAWGNQPGAGGYPGAAYP GAYPGQAPPGAYPGQAPPGA YPGQAPPSAYPGPTAPGAYP GPTAPGAYPGQPAPGAFFPGQ PGAPGAYPQCSGGYPAAGPY GVPAGPLTVPYDLPLPGGVM PRMLITIMGTVKPNANRIVL DFRRGNDVAFHFNPRFNENN RRVIVCNTKQDNNWGKEERQ SAFFFESGKPFKIQVLVEAD HFKVAVNDAHLLQYNHRMK
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.