

Mouse Amphiregulin/AREG protein (Recombinant) (N-hFc) (STJP020371)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11839
Gene Symbol	Areg
UniProt ID	AREG_MOUSE
Immunogen Region	Val100-Ala248
Immunogen Sequence	VIKPKKNKTEGEKSTKPKR KKKGGKNGKGRRNKKKKNPC TAKFQNFCHGECRYIENLE VVTCHNQDYFGERCGEKSM KTHSEDDKDLSKIYVAVTI FVSAILAAIGIGIVTVHL WKRYFREYEGETEERRRLRQ ENGTVHAIA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---