

Human disulfide-isomerase/PDI protein (Recombinant) (C-His) (STJP020372)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5034
Gene Symbol	P4HB
UniProt ID	PDIA1_HUMAN
Immunogen Region	Asp18-Lys505
Immunogen Sequence	DAPEEEDHVLVRKSNFAEA LAAHKYLLVEFYAPWCGHCK ALAPEYAKAAGKLKAEGSEI RLAKVDATEESDLAQQYGVR GYPTIKFFRNGDTASPKEYT AGREADDIVNWLKKRTGPAA TTLPDGAAAESLVESSEVAV IGFFKDVESDSAKQFLQAAE AIDDIPFGITSNSDVFSKYQ LDKDGVVLFKKFDEGRNNFE GEVTKENLLDFIKHNQLPLV IEFTEQTAPKIFGGEIKTH

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---