

## Human HGF Activator/HGFA protein (Recombinant) (C-6His) (STJP020374)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3083</a>
Gene Symbol	<a href="#">HGFA</a>
UniProt ID	<a href="#">HGFA_HUMAN</a>
Immunogen Region	Gln36-Ser655
Immunogen Sequence	QPGGNRTESPEPNATATPAI PTILVTSVTSETPATSAP EA EGPQSGGLPPPPRAVPSSSS PQAQALTEDGRPCRFPFRYG GRMLHACTSEGS AHRKWCAT THNYDRDRAWGYCVEATPPP GGPAALDPCASGPCLNGGSC SNTQDPQSYHCSPRAFTGK DCGTEKCFDETRYEYLEGGD RWARVRQGHVEQCECFGGRT WCEGTRHTACLSSPCLNGGT CHLIVATGTTVCACPPGFA

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---