

Human NUDT5 protein (Recombinant) (N-His) (STJP020376)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11164
Gene Symbol	NUDT5
UniProt ID	NUDT5_HUMAN
Immunogen Region	Met1-Phe219
Immunogen Sequence	MESQEPTESSQNGKQYIISE ELISEGKWVKLEKTTYMDPT GKTRTWESVKRTTRKEQTAD GVAVIPVLQRTLHYECIVLV KQFRPPMGGYCIEFFAGLID DGETPEAAALRELEEETGYK GDIAECSPAVCMDPGLSNCT IHIVTVTINGDDAENARPKP KPGDGEFVEVISLPKNDLLQ RLDALVAEEHLTVDARVYSY ALALKHANAKPFEVPFLKF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---