

## Human Histone H2A protein (Recombinant) (No tag) (STJP020377)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">92815</a>
Gene Symbol	<a href="#">H2AC25</a>
UniProt ID	<a href="#">H2A3_HUMAN</a>
Immunogen Region	Ser2-Lys130
Immunogen Sequence	SGRGKQGGKARAKAKSRSSR AGLQFPVGRVHLLRKGNYE ERVAGAPVYLAAVLEYLTA EILELAGNAARDNKKTRIIP RHLQLAIRNDEELNKLGRV TIAQGGVLPNIQAVLLPKKT ESHHKAKGK

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---