

Human Annexin A1/ANXA1 protein (Recombinant) (C-His) (STJP020378)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	301
Gene Symbol	ANXA1
UniProt ID	ANXA1_HUMAN
Immunogen Region	Met1-Asn346
Immunogen Sequence	MAMVSEFLKQAWFIENEEQE YVQTVKSSKGGPGSAVSPYP TFNPSSDVAALHKAIMVKGV DEATIIDILTKRNNAQRQQI KAAYLQETGKPLDETLKKAL TGHLEEVVLALLKTPAQFDA DELRAAMKGLGTDEDTLIEI LASRTNKEIRDINRVYREEL KRDIAKDITSDTSGDFRNAL LSLAKGDRSEDFGVNEDLAD SDARALYEAGERRKGTDVNV FNTILTTRSYPLRRVFQK
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.