

Human Inhibin beta E/INHBE protein (Recombinant) (N-hFC) (STJP020380)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83729
Gene Symbol	INHBE
UniProt ID	INHBE_HUMAN
Immunogen Region	Thr237-Ser350
Immunogen Sequence	TPTCEPATPLCCRRDHYVDF QELGWRDWILQPEGYQLNYC SGQCPPHLAGSPGIAASFHS AVFSLKANNWPASTSCCV PTARRPLSLLYLDHNGNVVK TDVPDMVVEACGCS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---