

Human IL-1 alpha protein (Recombinant) (NO-Tag) (STJP020382)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3552
Gene Symbol	IL1A
UniProt ID	IL1A_HUMAN
Immunogen Region	Ser113-Ala271
Immunogen Sequence	SAPFSFLSNVKYNFMRIIKY EFILNDALNQSII RANDQYL TAAALHNLDEAVKFDMGAYK SSKDDAKITVILRISKTLQY VTAQDEDQPVLLKEMPEIPK TITGSETNLLFFWETHGTKN YFTSVAHPNLFIA TKQDYWV CLAGGPPSITDFQILENQA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---