

Human IL-21 protein (Recombinant) (N-HIS & AVI) {Biotin} (STJP020383)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in 20mM NaAc, 500mM NaCl (pH 5.5). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	59067
Gene Symbol	IL21
UniProt ID	IL21_HUMAN
Immunogen Region	Gln32-Ser162
Immunogen Sequence	QDRHMIRMRLIDIVDQLKN YVNDLVPEFLPAPEDVETNC EWSAFSCFQKAQLKSANTGN NERIINVSIIKLLKRKPPSTN AGRRQKHRLTCPCDSYEKK PPKEFLERFKSLQKMIHQH LSSRTHGSEDS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---