

Human SLC3A2/CD98 protein (Recombinant) (N-His) (STJP020385)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	6520
Gene Symbol	SLC3A2
UniProt ID	4F2_HUMAN
Immunogen Region	Arg206-Ala630
Immunogen Sequence	RAPRCRELPAQKWWHTGALY RIGDLQAFQGHGAGNLAGLK GRLDYLSLKVKGVLVLP KNQKDDVAQTDLLQIDPNFG SKEDFDSLLQSAKKKSIRVI LDLTPTYRGENSWFSTQVDT VATKVKDALEFWLQAGVDGF QVRDIENLKDASSFLAEWQN ITKGFSEDRLLIAGTNSSDL QQILSLLESNKDLLTSSYL SDSGSTGEHTKSLVTQYLN TGNRWCSWSLSQARLLTSF
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.