

Mouse Anti-HIV Hiv antibody {Unconjugated} (STJS000002)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Mouse
Reactivity	HIV

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	WB 1:1000-10000 ELISA 1:10000-20000 IHC 1:50-100 IF 1:50-1000
Formulation	Liquid in PBS containing 0.5% BSA and 0.02% Sodium Azide.
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---