

Rabbit Anti-BSA BSA antibody {Biotin} (STJS000075)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	BSA

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Biotin
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	ELISA 1:10000-20000 IHC 1:50-200 IF 1:50-1000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---