

## Rabbit Anti-HSA HSA antibody {Biotin} (STJS000078)

### GENERAL INFORMATION

|               |                      |
|---------------|----------------------|
| Product Type  | Secondary antibodies |
| Applications  | IHC/IF/ELISA         |
| Host / Source | Rabbit               |
| Reactivity    | HSA                  |

### PRODUCT PROPERTIES

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                      |
| Concentration       | 2 mg/mL                                                                                         |
| Conjugation         | Biotin                                                                                          |
| Purification        | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range      | ELISA 1:10000-20000<br>IHC 1:50-200<br>IF 1:50-1000                                             |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                         |
| Isotype             | IgG                                                                                             |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.  |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|