

## Rabbit Anti-Donkey IgG Heavy & Light Chain antibody {HRP} (STJS000126)

### GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	ELISA/WB/IHC/ICC
Host / Source	Rabbit
Reactivity	Donkey

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	HRP
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	ELISA 1:5000-1:20000 WB 1:5000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---