

Goat Anti-Human IgM Mu Chain antibody {FITC} (STJS000181)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	IF/FC
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	FITC
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IF 1:250-1:2000 FCM 1:250-1:2000
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---