

Rabbit Anti-Guinea Pig IgG Heavy & Light Chain antibody {FITC} (STJS000193)

GENERAL INFORMATION

| | |
|---------------|----------------------|
| Product Type | Secondary antibodies |
| Applications | IF/FC |
| Host / Source | Rabbit |
| Reactivity | Guinea Pig |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Concentration | 1 mg/mL |
| Conjugation | FITC |
| Purification | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Dilution Range | IF 1:250-1:2000 FCM 1:250-1:2000 |
| Formulation | Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|