

Goat Anti-Mouse Mouse antibody {FITC} (STJS000794)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	IHC/IF
Host / Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	FITC
Purification	IgG fraction purified with Protein A affinity chromatography
Dilution Range	IF: 1 to 2 μ g/ml Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 50% glycerole
Isotype	IgG
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Specificity	Reacts with mouse IgG-s
-------------	-------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---