

Goat Anti-Rat IgG (H&L) antibody {Alexa Fluor 488} (STJS000814)

GENERAL INFORMATION

Product Type	Secondary antibodies
Host / Source	Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Alexa Fluor 488
Purification	Affinity purified
Dilution Range	Fluorescence microscopy: 1-2 ug/mL, Flow cytometry: 1 ug/106 cells, Near-infrared western detection: 50-100 ng/mL
Formulation	Yellow solution containing pH~7.4 PBS containing 50 mg/mL glycerol 2 mg/mL bovine serum albumin (IgG-free and protease-free) and 0.0.5% sodium azide.
Isotype	IgG
Storage Instruction	Store at 4°C for several months. Protect from light. For longer storage, divide the conjugate into small aliquots and freeze at -20°C. Avoid repeated freezing and thawing.

TARGET INFORMATION

Immunogen Region	H&L
------------------	-----

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---