

Goat Anti-Mouse IgG Highly Cross-Adsorbed (H&L) antibody {Alexa Fluor 647} (STJS000848)**GENERAL INFORMATION**

Product Type	Secondary antibodies
Host / Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Alexa Fluor 647
Purification	Affinity purified
Dilution Range	Fluorescence microscopy: 1-2 ug/mL. Flow cytometry: 1 ug/106 cells. Near-infrared western detection: 50-100 ng/mL
Formulation	Blue solution in pH~7.4 PBS containing 50% glycerol, 2 mg/ml bovine serum albumin (IgG-free and protease-free) and 0.05% sodium azide.
Isotype	IgG
Storage Instruction	Store at 4°C for several months. Protect from light. For longer storage, divide the conjugate into small aliquots and freeze at -20°C. Avoid repeated freezing and thawing.

TARGET INFORMATION

Immunogen Region	H&L
------------------	-----

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---