

Goat Anti-Rat IgG H&L antibody {HRP} (STJS000857)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/IHC-P
Host / Source	Goat
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB:1:5000-1:10000 IHC-P:1:500-1:2000
Formulation	PBS with 0.75% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	Rat IgG
-------------	---------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of Rat IgG, using HRP Goat Anti-Rat IgG (H+L) antibody (STJS000857) at 1:1000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.