

Mouse Anti-Duck IgY H&L antibody {AMCA} (STJS000894)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	ELISA/IHC/IF/FC
Host / Source	Mouse
Reactivity	Duck

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	AMCA
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgY
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---