

**Goat Anti-Chicken IgY H&L antibody {TRITC} (STJS000906)****GENERAL INFORMATION**

|               |                      |
|---------------|----------------------|
| Product Type  | Secondary antibodies |
| Applications  | ELISA/IHC/IF/FC      |
| Host / Source | Goat                 |
| Reactivity    | Chicken              |

**PRODUCT PROPERTIES**

|                     |                                                                                                                   |
|---------------------|-------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                        |
| Conjugation         | TRITC                                                                                                             |
| Purification        | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Dilution Range      | IHC 1:200-1:1000<br>IF 1:200-1:1000<br>FCM 1:2000-1:4000<br>ELISA                                                 |
| Formulation         | Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide                                            |
| Isotype             | IgY                                                                                                               |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                    |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|