

Rabbit Anti-Sheep IgG H+L antibody {HRP} (STJS001174)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Sheep

PRODUCT PROPERTIES

Concentration	Lot specific
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB:1:5000-1:100000 IHC-P:1:500-1:5000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.75% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	Sheep IgG
-------------	-----------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---