

Rabbit Anti-Goat IgG H+L antibody {HRP} (STJS001178)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/IHC-P
Host / Source	Rabbit
Reactivity	Goat

PRODUCT PROPERTIES

Concentration	Lot specific
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB:1:5000-1:100000 IHC-P:1:500-1:5000
Formulation	PBS with 0.025% Sodium Azide, 0.75% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	Goat IgG
-------------	----------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---