

Goat Anti-Rabbit IgG H+L antibody {ABflo® 488} (STJS001206)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	IF/ICC/FC
Host / Source	Goat
Reactivity	Rabbit

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	ABflo® 488
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 FC:1:100-1:800
Formulation	PBS with 0.025% Sodium Azide, 0.75% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Specificity	Rabbit IgG
-------------	------------

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




Immunofluorescence analysis of HeLa cells using Alpha-tubulin (ubiquitous) chain Rabbit mAb (STJ11103694) (40x lens), the secondary antibody was ABflo® 488-conjugated AffiniPure Goat Anti antibody (STJS001206) used at dilution of 1 : 200. Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (STJS001166) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using Alpha-tubulin (ubiquitous) chain Rabbit mAb (STJ11103694) (40x lens), the secondary antibody was ABflo® 488-conjugated AffiniPure Goat Anti antibody (STJS001206) used at dilution of 1 : 200. Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (STJS001166) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of PC-12 cells using Alpha-tubulin (ubiquitous) chain Rabbit mAb (STJ11103694) (40x lens), the secondary antibody was ABflo® 488-conjugated AffiniPure Goat Anti antibody (STJS001206) used at dilution of 1 : 200. Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (STJS001166) at 1:500 dilution. Blue: DAPI for nuclear staining.



Flow cytometric analysis of Positive antibody E-Cadherin Rabbit mAb (2.5 µg/mL) in various cells (orange) compare to Rabbit rabbit isotype control (blue) and non-staining control (Red). The secondary antibody used was ABflo® 488-conjugated AffiniPure Goat Anti-Rabbit IgG (H+L) (STJS001206) at 1:100.