

Rabbit Anti-Human IgG1 Fc antibody {HRP} (STJS001216)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB:1:2000-1:30000 ELISA: 1:5000-1:10000
Formulation	PBS with 0.03% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	50kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	Human IgG
-------------	-----------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates, using HRP Rabbit anti-Human IgG1 (Fc) (STJS001216) at 1:4000-1:32000 dilution.

Lysates/proteins: 25 µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit

Exposure time: 30s.