

Mouse Anti-Rabbit Rabbit IgG antibody {HRP} (STJS001235)

GENERAL INFORMATION


Product Type	Secondary antibodies
Applications	WB/IP/ELISA
Host / Source	Mouse
Reactivity	Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity Purified
Dilution Range	WB-1:2500-1:10000 ELISA-1:1000-1:3000
Formulation	Liquid in phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, 50% glycerol, pH 7.4.
Isotype	IgG2b
Storage Instruction	Protect from light and store at 4°C short term. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti-Rabbit IgG antibody western blot native
 denatured rabbit IgG STJS001235