

Rat Anti-Mouse Mouse IgG antibody {HRP} (STJS001236)

GENERAL INFORMATION


Product Type	Secondary antibodies
Applications	WB/IP/ELISA
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity Purified
Dilution Range	WB-1:2500-1:10000 ELISA-1:1000-1:3000
Formulation	Liquid in phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, 50% glycerol, pH 7.4. NA/mL
Storage Instruction	Protect from light and store at 4°C short term. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IgG antibody western blot mouse IgG native
denatured HRP STJS001236