

**Mouse Anti-Human Mouse Human IgG Lambda light chain antibody {Unconjugated}
(STJS001237)**

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/IHC-F/IHC-P/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Formulation	Liquid in 0.01M PBS, PH:7.2, Preservative.
Storage Instruction	Store at-20°C and avoid freeze/thaw cycles. Avoid light.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---