

Mouse Anti-Species-independent His Tag antibody {HRP} (STJS001238)

GENERAL INFORMATION


Product Type	Secondary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Species-independent

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity Purified
Dilution Range	WB 1:1000-1:10000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



His Tag antibody western blot recombinant protein STJS001238



His Tag antibody western blot 12tag protein STJS001238