

Mouse Anti-Species-independent Flag Tag antibody {HRP} (STJS001239)

GENERAL INFORMATION


Product Type	Secondary antibodies
Applications	WB/IP
Host / Source	Mouse
Reactivity	Species-independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	HRP
Purification	Affinity Purified
Dilution Range	WB 1:1000-1:10000 IP 1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 2-8°C. Do not freeze.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Flag tag antibody western blot CHO-K1 cells
 STJS001239