

Mouse Anti-Species-independent Flag Tag antibody {HRP} (STJS001241)

GENERAL INFORMATION


Product Type	Secondary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Species-independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity Purified
Dilution Range	WB 1:1000-1:10000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



FLAG tag antibody western blot CHO-K1 cells
 STJS001241