

Goat Anti-Chicken IgY antibody {FITC} (STJS001244)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	IHC/ICC/WB
Host / Source	Goat
Reactivity	Chicken

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	Phosphate-buffered saline (PBS) (pH 7.2) , bovine serum albumin (0.2%, w/v) , sodium azide (0.02%, w/v) (as an anti-microbial agent)
Isotype	IgG
Storage Instruction	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. For longer term storage, aliquot and freeze to avoid freeze-thaw of the antibody.

TARGET INFORMATION


Immunogen Sequence	Chicken
Specificity	H&L; does not cross-react with mammalian IgGs.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



IHC staining rat brain cerebellum using goat anti-chicken IgY FITC secondary antibody STJS001244



IHC staining rat brain cerebellum using goat anti-chicken IgY FITC secondary antibody STJS001244