

Chicken Anti-Rabbit IgG antibody {FITC} (STJS001251)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	IHC/ICC/WB
Host / Source	Chicken
Reactivity	Rabbit

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	FITC
Purification	Affinity Purified
Formulation	Phosphate-buffered saline (PBS) (pH 7.2) , sodium azide (0.02%, w/v) (as an anti-microbial agent)
Isotype	IgY
Storage Instruction	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. For longer term storage, aliquot and freeze to avoid freeze-thaw of the antibody.

TARGET INFORMATION

Immunogen Sequence	Rabbit
Specificity	H&L; there is no detectable cross-reactivity with goat, human or mouse IgG.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---