

Chicken Anti-Rabbit IgG antibody {HRP} (STJS001252)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	IHC/ICC/WB
Host / Source	Chicken
Reactivity	Rabbit

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	HRP
Purification	Affinity Purified
Formulation	HEPES (50 mM, pH 7.0) , NaCl (0.1 M) , bovine serum albumin (0.2%, w/v) , thimerosal (0.01%, w/v) (as an anti-microbial agent)
Isotype	IgY
Storage Instruction	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. For longer term storage, aliquot and freeze to avoid freeze-thaw of the antibody.

TARGET INFORMATION

Immunogen Sequence	Rabbit
Specificity	H&L; there is no detectable cross-reactivity with goat, human or mouse IgG.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---