

Rabbit Anti-Human IgA Fc antibody {HRP} (STJS001264)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity purification
Dilution Range	ELISA, 1:5000-1:10000
Isotype	IgG
Molecular Weight	Protein Weight Observed Mw 83kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	IGHA1 Uniprot ID= IGHA1_HUMAN "
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---